

# Nox A1 - Special Features and Benefits

## Customer Support Document

Version 1.0

Latest Revision 4 November 2013

Copyright© 2013

Nox Medical - All rights reserved

Nox A1<sup>®</sup>, Nox T3<sup>®</sup> and Noxturnal<sup>®</sup> are manufactured by:

Nox Medical ehf

Katrinartuni 2

IS - 105 Reykjavik

Iceland

Website: [www.noxmedical.com](http://www.noxmedical.com)



### Copyright Notice

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language or computer language, in any form, or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without the prior written authorization from Nox Medical.

### Disclaimer

This document may contain typographical errors or technical inaccuracies. Nox Medical does not accept any liability for the use or misuse whether direct or indirect of the products, or for damages arising out of the use of or inability to use the products. Users must accept all responsibility for any results obtained by or concluded from data obtained by the products including software from Nox Medical. All clinical conclusions and decisions that are based on the use of this product are the responsibility of the user.

## Overview

This knowledge base article describes the special features and benefits of the Nox A1 PSG System. It is intended to help users understand the true value of using the A1 system and how it can reflect on their daily work.

## Common Issues Being Solved With the Nox A1 PSG System

The company's mission is to improve health by improving sleep diagnostic by making it more simple, efficient and comfortable in all patient groups, especially children.

In our aim to simplify sleep diagnostics, our engineering team identified multiple issues with traditional PSG systems which are available in the field. The intention of releasing the A1 system was to address many of these issues to advance sleep medicine in general. Examples of common issues with traditional PSG systems are:

- **Equipment size and bulkiness**  
This makes many PSG systems uncomfortable for patients and impractical for pediatric use.
- **Tethered patients**  
In many cases, patients need to be tethered to the bed with electrodes, wires, and cables. This makes many patients uncomfortable and makes it difficult for them to use the restroom at night.
- **Cables get tangled**  
Many PSG systems rely on long cables which easily get in the way of patients natural flow of movement. This makes patients uncomfortable and increases the risk of cables getting tangled and even pulled off the patients, increasing risk of lost studies.
- **Lost studies**  
Do to complexities in setting up many traditional PSG systems and their reliance on a few key sensors, many studies are lost. Electrodes may fall off, the signal quality may deteriorate throughout the night or patients might lose sleep because of discomfort. Redoing sleep studies can be very expensive.

The Nox A1 PSG system addresses all of these common issues with multiple features. Below you will find a list of the key features and their benefits. After reviewing this list you will be able to understand the key benefits of the Nox A1 PSG System.

## NOX PSG – special features and benefits

Feature	Benefit	Target
Smallest and most lightweight PSG recorder	Patient comfort and flexibility for ambulatory and pediatric studies → study patients in their natural sleep environment	Patient comfort, System flexibility
Reduced complexity - simplified PSG setup	Time and effort saving	Reduced hookup time
Ergonomic setup with minimum number of cables and straps	Minimum of patient impact → comfort for study acceptance, easy setup enabling home sleep testing	Patient comfort
Unique EEG leads setup	One multi-lead cable from recorder <ul style="list-style-type: none"> <li>- minimize electrode pull → reduced risk of lost EEG electrodes</li> <li>- minimize cable tangling</li> <li>- intuitive lead choice / placement</li> </ul>	Study integrity, Patient comfort, Reduced hookup time
New ECG and EMG leg movement cables	Double lead structure to minimize number of single cables and quicker setup	Patient comfort, Reduced hookup time
EEG and chin EMG: <ul style="list-style-type: none"> <li>• Continuous impedance control</li> <li>• Independent floating reference (automatic reference adjustment)</li> </ul>	High signal quality for all well conducted EEG / EMG signals independent from any impedance issues or lost electrodes	Signal quality, Study integrity
Simultaneous digitalization for all channels	Signal synchrony, minimizing crosstalk and noise introduction by sampling skewing	Signal quality
Sampling frequency 256 kHz	Extreme oversampling for high-end noise reduction and anti-aliasing, allow clean signals free of signal disturbance	Signal quality
Floating PSG system	No power cable → no power-line disturbance → no need for extra filtering and lost information e.g. in EEG based on bandpass stop filters	Signal quality → especially complete EEG signal information
32 bit signal processing (16 bit dynamic storage)	High signal quality by minimizing quantification noise and by using arithmetic in highest available resolution before storing	Signal quality
Build-in calibrated Microphone	Recording of calibrated sound – dbA (decibel audio) similar to high end microphones – allows true sound playback and detailed snoring frequency analysis	Reliable snoring analysis
High quality RIP allowing calibrated RIP technology	Current RIP technology usually used for backup flow calculation based on correction for phase angle, calibrated RIP would replace such correction by a true flow calculation (also offering flow volume loop visualization) → enhanced flow analysis – esp. UARS, enhanced titration options, reliable backup flow	Signal quality for enhanced analysis (diagnostic and titration), Study integrity
Integrated ambient light sensor	Light off / on detection for easier analysis, especially for ambulatory studies	Scoring support
Bipolar channels with programmable current source	Possibility to switch from passive bipolar sensor recording to actively drive special sensors (Galvanic Skin Response Sensors – night sweating recording)	Flexible system use

## Hardware related software features

Feature	Benefit	Target
Nox calibrated RIP Flow (based on high quality RIP technology)	Reliable flow calculation to analyze the complete respiratory system breathing information (independent of nose or mouth breather), enhanced apnea, hypopnea and flow limitation scoring	Reliable Flow Analysis, Study integrity
Flow Volume Loops (based on calibrated RIP Flow derived on high quality RIP technology)	Improved analysis of UARS, extended titration possibilities to titrate back to a signature breath	Enhanced analysis, extended parameters during titration
Video integration (synchronized) for ambulatory studies (based on specific available signal information)	Easy setup video for any kind of ambulatory study setup (HST or hospital)	Flexibility on HST
Tablet application (based on build in Bluetooth)	Direct study setup and signal control	Flexibility on HST, study success rate

## Special PSG software features

Feature	Benefit	Target
Unique PSG related "Recording Results Page"	Easy complete PSG data overview for physician and technician	Reduced time for quick study overview
Extremely powerful workspace and scoring handling	Prepared customized views to load / change quickly between workspaces, split screens, scoring screens and overview windows	Flexibility for scoring and analysis views
Powerful and flexible reporting engine	Design reports like Word documents, easy calculation of new parameters from existing variables, reporting of various customizable event and/or signal associations, load existing Word reports into Noxturnal	Flexibility for reporting
Automatic event control	Avoid scoring or reporting on concurrent events	Accuracy on scoring and reporting
Automatic sleep staging	Allow quick sleep staging by automatic pre-staging and manual correction	Reduced time for quick analysis
Simplified user and administrator settings	Software setup customization and protection, easy portability of custom settings	Reduced installation time, enhanced user customization and settings control
Data export as signal / event grids	Export structures can customized and saved as sheet templates	Scientific use: export signal and event information