

Designed to easily fit them and everyone in between



AirFit™ 20 series

Thanks to the InfinitySeal™ cushion,
AirFit F20 fits 96,5% and
AirFit N20, 99,4% of patients tested.¹⁻⁵
They come in 'for Her' variants too!
All are made to seal robustly
regardless of facial shape or size.

- 1 ResMed internal study of 22 existing ResMed patients, conducted between 26/04/2016 – 27/05/2016 comparing the current market leading mask with AirFit F20. Preliminary patient study – data on file, ID A3810791.
- 2 ResMed AirFit F20 internal fitting study of 27 ResMed and non-ResMed patients, conducted between 30/03/2016 – 04/04/2016. Preliminary patient study – data on file, ID A3751086.
- 3 ResMed AirFit F20 internal international fitting study of 34 ResMed and non-ResMed patients, conducted between 11/04/2016 – 15/04/2016. Preliminary patient study – data on file, ID A3774922.
- 4 ResMed AirFit F20 internal international fitting study of 90 ResMed and non-ResMed patients conducted between 06/06/2016 – 22/06/2016. Preliminary patient study – data on file, ID A3430701.
- 5 ResMed AirFit N20 internal international fitting study of 159 existing Resmed patients, conducted 12/11/2015; ID A3697629.

Find out more at ResMed.com/AirFit20

© 2016 ResMed Ltd. 102259/1 2016-11



ResMed